

XF-1600™ Flow Cytometer

Expert Technology
for Laboratories



Together for a better
healthcare journey

For Research Use Only. Not for use in diagnostic procedures.
RUO instruments must be validated before use in clinical practice.

XF-1600 Flow Cytometer

Sysmex, a global leader in hematology, urinalysis and hemostasis and a provider of trusted automated workflow solutions, now offers scientists a unique and innovative flow cytometry solution for robust data collection.

The XF-1600 integrates a multi-laser optical layout with Sysmex's proven fluidics design for a reliable flow cytometer performance. Stable fluidics—even at high sample acquisition rates—ensure the system is capable of rapid data collection and analysis with high sensitivity. The XF-1600 also provides the option of a lower sampling rate for increased measurement precision.

Designed for increased stability



Flexible specimen handling

Rare event analysis capability

Rapid data acquisition without loss of sensitivity

Intuitive software with automated QC and optional off-line analysis solution

Top-rated service sysmex is known for

Technology meets reliability



The XF-1600 uses innovative technology and fluidics to provide high detection sensitivity, even at high sheath velocity, allowing for rapid data acquisition. Like Sysmex's clinical testing platforms, the XF-1600 flow cytometer is backed by the exceptional service and reliability laboratories expect from a global leader in hematology.

- 37 measurable parameters simultaneously record signal-area, -height, -width and -time data from 12 detection channels
- Linear, log and logicle plots are available for data acquisition
- Complete array of quality control products with automated QC setup
- Proven Sysmex fluidics design

For Research Use Only. Not for use in diagnostic procedures. RUO instruments must be validated before use in clinical practice.

Powerful Rare Event Mode with excellent data resolution

With the growing interest for rare event analysis, the XF-1600 incorporates a Rare Event Mode (REM) based on established hematology technology. The system allows users a choice of six aspiration volumes which analyze up to 2 mL of specimen for low carryover and Limit of Blank (LoB).




Five adjustable sample rates - from 5 $\mu\text{L}/\text{min}$ to 120 $\mu\text{L}/\text{min}$ - provide flexibility for all sample types, titer levels and cell concentrations. At high sheath velocity, fast data acquisition rates - up to 50,000 events per second - reduce the overall sample acquisition time while maintaining sensitivity required for rare event analysis.

Data acquisition rates up to 50,000 events per second reduce overall analysis time and provide an impressive 5 million events of robust data.

8 Peak Beads at 5 $\mu\text{L}/\text{minute}$ flow rate




8 Peak Beads at 120 $\mu\text{L}/\text{minute}$ flow rate



For Research Use Only. Not for use in diagnostic procedures. RUO instruments must be validated before use in clinical practice.


Intuitive XF-1600 software

Sysmex's intuitive user interface makes the XF-1600 easy to use for all lab staff. Choose from a variety of plot options including linear, log and logicle axis scaling for single or dual parameter plot types with selectable statistics.



Data is for representative purposes only. Complex tests must be validated by laboratories.

For Research Use Only. Not for use in diagnostic procedures. RUO instruments must be validated before use in clinical practice.



Variety of options for panel customization:

- Easy selection and input of channel settings
- Five different flow rates
- Volume selection of up to 2 mL
- Multiple stop settings based off event count or time

Channel Setting

LOGIC ☐ OR ☒ AND

Channel	Name	Height	Area	Width	Threshold	Gain/PMT
<input checked="" type="checkbox"/>	FSC	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	2500	333
<input checked="" type="checkbox"/>	SSC	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	1	515
<input type="checkbox"/>	FL1 CD2-FITC	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	628
<input type="checkbox"/>	FL2 CD22-PE	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	585
<input type="checkbox"/>	FL3 CD16-PEDL	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	655
<input type="checkbox"/>	FL4 CD19-PERCP	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	665
<input type="checkbox"/>	FL5 CD8-PECY7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	658
<input type="checkbox"/>	FL6 CD56-APC	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	785
<input type="checkbox"/>	FL7 CD3-AF700	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	730
<input type="checkbox"/>	FL8 CD4-APCY7	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	680
<input type="checkbox"/>	FL9 CD20-PACB	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	504
<input type="checkbox"/>	FL10 CD45-PACO	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	1	575

Flow Rate ☐ High ☒ Middle ☐ Low
 60


Volume Required volume select 500 uL

Laser 488 638 405 [nm]

Stop Setting
 Event Count ☒ Stop Condition Ungated Event Count 30000
 Time ☐ Stop Time 68 Seconds (5s-226s)

For Research Use Only. Not for use in diagnostic procedures. RUO instruments must be validated before use in clinical practice.

A flexible approach to sample handling



XF-1600 shown with Helmer 24-tube rotor

If your valuable specimens require special attention, the XF-1600 provides unique sampling options. Single tube processing within the manual tube position can be performed on a variety of tube sizes with selectable aspiration volumes to ensure minimal sample waste and protection against accidental sample loss. Alternatively, samples can be loaded on a Sysmex compatible rotor for automated walkaway acquisition.

For Research Use Only. Not for use in diagnostic procedures. RUO instruments must be validated before use in clinical practice.



Seamless integration

For exceptional efficiency, Sysmex offers additional options for streamlining workflow in the flow cytometry lab.

PS-10™ Sample Preparation System

The PS-10 Sample Preparation System automates many manual tasks associated with preparing flow cytometry samples, freeing up valuable time for techs. Additionally, it instills trust and confidence in data due to its excellent sample viability, high cell recovery and reproducible results.

Helmer UltraCW® II Automatic Cell Washing System

For procedures that include wash steps, the Helmer UltraCW II Automatic Cell Washing System integrates seamlessly into your flow lab processes. The barcoded rotors are compatible with the XF-1600 and PS-10, offering positive identification to ensure traceability through the laboratory.

VenturiOne® Software




VenturiOne offline analysis software provides high-speed processing of up to 400 data files, with unique and robust data previews that make plot creation simple and quick. Autogating, hyperlog slides and versatile color compensation put powerful analytics in your hands.

Specifications

Principles and Technology Standard single tube loading port and rotor-based autoloader with optional sample mixing capabilities
3 spatially separated lasers: 405 nm, 488 nm, and 638 nm
Operating system (OS): Windows® 10 professional (64bit)
FCS format 3.0

Fluorescence Detectors 10 fluorescent channels

Laser	Channel	Filter	Laser	Channel	Filter
488 nm	FL1	530/30	638 nm	FL6	660/20
	FL2	580/30		FL7	720/40
	FL3	610/14		FL8	780/60
	FL4	700/50	405 nm	FL9	450/50
	FL5	750 LP		FL10	540/80



Sample Flow Rate
Low: 5 µL per minute
Medium: 30 µL per minute and 60 µL per minute
High: 90 µL per minute and 120 µL per minute

Sensitivity: MESF
FITC < 100
PE < 50
APC < 100
Particle Resolution - 0.5 µm on scatter and up to 50 µm

Features
Fast data acquisition rates up to 50,000 events per second and 5 million events from 37 parameters including height, area, width and time
Carry-over at 0.1% in Standard mode

PS-10 compatible

Dimensions and Weight
Approx. 27.5 (W) X 25.0 (D) X 29.0 (H) in., approx. 220 lbs
700 (W) x 630 (D) x 740 (H) mm, approx. 100 kg



Award Winning Service, Reliability and Training

For more than half a century, Sysmex has been providing diagnostic equipment to clinical labs around the globe. Years of experience utilizing flow cytometry principles in hematology and urinalysis instrumentation have led to the development of the XF-1600, our first flow cytometer for research use.

Sysmex is well known for providing best-in-class service and reliability alongside our superior technology. A robust training curriculum, including e-learning and live virtual instructor-led training (VILT), is provided by the award-winning Sysmex Center for Learning. Additionally, Sysmex's hematology systems consistently earn the Best Service, Best Customer Satisfaction and Best System Performance awards* from IMV ServiceTrak™ and nine of the top ten cancer care hospitals in the United States use Sysmex analyzers.†

*Source: <https://imvinfo.com/press-room/imv-announces-the-first-wave-of-2023-imv-servicetrak-clinical-laboratory-awardees-at-aacc-2023/>

†Source: <https://health.usnews.com/best-hospitals>

For the Technical Assistance Center (TAC), call 1-888-879-7639. For Canada TAC, call 1-888-679-7639.

For Research Use Only. Not for use in diagnostic procedures. RUO instruments must be validated before use in clinical practice.

Sysmex Corporation
1-5-1-Wakoinohama-Kaigandori Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0521 · www.sysmex.co.jp

Sysmex America, Inc.
577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A. · Phone +1 800 379-7639 · www.sysmex.com/us

Sysmex Canada, Inc.
5700 Explorer Drive Suite 200, Mississauga, ON L4W0C6 Canada · Phone +1 905 366-7900 · www.sysmex.ca

Sysmex Latin America and the Caribbean
Rua Joaquim Nabuco 615 - Bairro Cidade Jardim, São José dos Pinhais Paraná - Brasil - CEP 83040-210 · Phone +55 41 2104-1314 · www.sysmex.com.br