

A CLIA-WAIVED CBC  
IS NOW  
POSSIBLE.



CLIA Waived  
XW-100



RxOnly





## TECHNOLOGY IS EVOLVING

Tests that once required a central lab can now be performed at the point-of-care. These changes are benefitting both patients and providers through faster results which can:

- Expedite diagnosis and treatment
- Improve patient satisfaction
- Streamline workflow

## A HISTORY OF FIRSTS

We are the global leader in hematology analyzers having introduced our 1st blood counter in 1963 and the first automated hematology system in 1988.

Today, Sysmex continues to meet the high standards that healthcare professionals require and expect.

## INTRODUCING THE SYSMEX® XW™-100

Our latest innovation is the Sysmex XW-100. It is the first FDA-cleared, CLIA-waived CBC analyzer to provide reliable, convenient, and often, same-visit CBC results. A 15 µL venous blood sample is required. The sample-to-result time is just 3 minutes.



## IT MAY FIT WELL IN YOUR PRACTICE

Previously, CBC testing required sending samples to a lab for results. The Sysmex XW-100 has changed that. The CLIA-waived designation ensures that it's simple to use, has a low risk of providing erroneous results, and can be operated without additional training beyond simply reading the manufacturer's instructions and following the on-screen prompts.

The Sysmex XW-100 can be an especially good fit for your well patient visits. It is very compact with a height of 13.8 inches and a width of 7.3 inches. The Sysmex XW-100 and its reagents can fit on a countertop. Daily QC takes less than 30 minutes.



## VALUABLE INFORMATION

The Sysmex XW-100 offers a 3-part differential with 12 different parameters:

- Total #WBCs
- Total #RBCs
- Hemoglobin
- Hematocrit
- Total #platelets
- Total #neutrophils
- % of neutrophils
- Total #lymphocytes
- % of lymphocytes
- Total #other WBCs
- % of other WBCs
- MCV

The Sysmex XW-100 is not for use in diagnosing or monitoring patients with primary or secondary chronic hematologic diseases/disorders, oncology patients, critically ill patients, or children under the age of two.

# SAME-VISIT CBC. IT IS POSSIBLE.


1 Prior to or during the examination, the patient's blood is drawn and the CBC is performed in your office.



2 CBC results are now available to discuss with your patient which allows for discussions and treatment decisions on the same visit.



3 Obtaining results with the XW-100 can help improve efficiency which may ultimately improve your bottom line. It may also increase patient satisfaction.



# PERFORMS EXACTLY AS DESIGNED

The Sysmex XW-100 reports the 12 parameters and normal reference ranges. It is a CLIA-waived analyzer that is factory calibrated and it will not report values that might be inaccurate.

SYSMEX XW-100 RESULTS (NON-SUPPRESSED)	
Instrument type XW100	
Serial # G2884	
Date Jan 14, 2019	
Time 10:21 AM	
Operator	DS
Patient ID	5832
Patient DOB	Aug 24, 1946
WBC	6.7 × 10 <sup>3</sup> /μL
RBC	4.44 × 10 <sup>6</sup> /μL
HGB	12.3 g/dL
HCT	37.1 %
PLT	227 × 10 <sup>3</sup> /μL
#Neut	3.9 × 10 <sup>3</sup> /μL
%Neut	58.4 %
#Lymph	2.2 × 10 <sup>3</sup> /μL
%Lymph	32.3 %
#OtherWBC	0.6 × 10 <sup>3</sup> /μL
%OtherWBC	9.3 %
MCV	83.6 fL
NOTES	
-----	
Adult Reference Ranges	
WBC	3.9 - 10.4 × 10 <sup>3</sup> /μL
RBC	3.71 - 5.52 × 10 <sup>6</sup> /μL
HGB	10.9 - 16.7 g/dL
HCT	32.5 - 49.4 %
PLT	148 - 382 × 10 <sup>3</sup> /μL
#Neut	2.2 - 7.1 × 10 <sup>3</sup> /μL
%Neut	46.4 - 76.9 %
#Lymph	0.9 - 3.4 × 10 <sup>3</sup> /μL
%Lymph	14.7 - 45.9 %
#OtherWBC	0.2 - 1.2 × 10 <sup>3</sup> /μL
%OtherWBC	3.2 - 16.9 %
MCV	82.5 - 98.0 fL
-----End Report-----	

PRACTICE EFFICIENCY +  
PATIENT SATISFACTION



## PROTECTING YOU AND YOUR PATIENTS

Blood parameters can be complicated to measure. The complexity of the sample and underlying patient conditions may result in suppression of results. This will appear as 4 asterisks (\*\*\*\*) where in most cases a result would be generated. The Sysmex XW-100 is designed to protect your patients and your practice from inaccurate results.

For more information, review the Sysmex XW-100 Quick Guide or visit [CBCin3.com](http://CBCin3.com).

**SYSMEX XW-100 RESULTS**  
 (SUPPRESSED)

Instrument type XW100  
Serial # G2883

Date Jan 14, 2019  
Time 12:16 PM  
Operator MKZ

Patient ID 1675  
Patient DOB May 26, 1972

WBC	6.2 × 10 <sup>3</sup> /μL
RBC	4.36 × 10 <sup>6</sup> /μL
HGB	****
HCT	****
PLT	344 × 10 <sup>3</sup> /μL
#Neut	4.4 × 10 <sup>3</sup> /μL
%Neut	71.3 %
#Lymph	1.6 × 10 <sup>3</sup> /μL
%Lymph	25.6 %
#OtherWBC	0.2 × 10 <sup>3</sup> /μL
%OtherWBC	3.1 % Low
MCV	****


**NOTES**

RECOMMEND FURTHER TESTING.

Adult Reference Ranges

WBC	3.9 - 10.4 × 10 <sup>3</sup> /μL
RBC	3.71 - 5.52 × 10 <sup>6</sup> /μL
HGB	10.9 - 16.7 g/dL
HCT	32.5 - 49.4 %
PLT	148 - 382 × 10 <sup>3</sup> /μL
#Neut	2.2 - 7.1 × 10 <sup>3</sup> /μL
%Neut	46.4 - 76.9 %
#Lymph	0.9 - 3.4 × 10 <sup>3</sup> /μL
%Lymph	14.7 - 45.9 %
#OtherWBC	0.2 - 1.2 × 10 <sup>3</sup> /μL
%OtherWBC	3.2 - 16.9 %
MCV	82.5 - 98.0 fL

-----End Report-----



Rerun sample if device alerts to do so. If results are still suppressed, send sample out as per your standard protocol.

## SYSMEX XW-100 INTENDED USE

The Sysmex XW-100 Automated Hematology Analyzer (XW-100) is a quantitative automated hematology analyzer intended for in vitro diagnostic use to classify and enumerate the following parameters for venous whole blood anticoagulated with K2/K3 EDTA: WBC, RBC, HGB, HCT, PLT, #NEUT, %NEUT, #LYM, %LYM, #OTHER WBC, %OTHER WBC, MCV. It is not for use in diagnosing or monitoring patients with primary or secondary chronic hematologic diseases/disorders, oncology patients, critically ill patients, or children under the age of two.

## SPECIFICATIONS

<b>Test Principle:</b>	DC (direct current) with hydrodynamic focusing for all parameters except hemoglobin, which is measured photometrically
<b>Sample Type:</b>	Venous whole blood
<b>Sample Volume:</b>	Approximately 15 μL (1mL of sample in the tube and 15μL to aspirate)
<b>Data Storage:</b>	100 test cycles
<b>Size and Weight:</b>	7.3" (W) x 18.1" (D) x 13.8" (H); 38 lbs
<b>Reagents:</b>	XW pack D (2 X 2L); XW pack L (2 X 250mL)
<b>Storage:</b>	All XW-100 reagents, XW CELLCLEAN™ and XW QC CHECK should be stored at room temperature (59° - 77° F or 15 - 25° C).
<b>Quality Control:</b>	Three levels of QC must be performed as follows: <ul style="list-style-type: none"> <li>• Every 8 hours of instrument use</li> <li>• When a new lot of reagent is used</li> <li>• After completing weekly care</li> </ul> The system will prompt users when QC is required and will not allow patient samples to be tested until QC is performed.

### Sysmex Part Numbers and Description

Material Number	Product
XW100-CALIB	XW-100 CALIBRATED
BL551785	XW CELLCLEAN (20 vials)
224032	XW QC CHECK (6X2.0mL)
ZA900005	XW pack D (2X2,000mL)
ZA900006	XW pack L (2X 250mL set)

### Consumables

Material Number	Product
HEMPOCHIPAPER	.25" x 55" white thermal paper 3 pack
44230318	SAMPLE TUBE ADAPTER (CREAM/13MM DIA)
44230512	SAMPLE TUBE ADAPTER (GREEN/CONTROL)
02324415	CUBITAINER SPOUT KIT NO.5 (fits D pack)
03304111	CUBITAINER SPOUT KIT NO.7 (fits L pack)

**CPT Code: 85025QW**



**Sysmex Corporation**  
1-5-1-Wakinohama-Kaigandori Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0521 · [www.sysmex.co.jp](http://www.sysmex.co.jp)

**Sysmex America, Inc.**  
577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A. · Phone +1 800 379-7639 · [www.sysmex.com/us](http://www.sysmex.com/us)

**Sysmex Canada, Inc.**  
5700 Explorer Drive Suite 200, Mississauga, ON L4W0C6 Canada · Phone +1 905 366-7900 · [www.sysmex.ca](http://www.sysmex.ca)

**Sysmex Latin America and the Caribbean**  
Rua Joaquim Nabuco 615 - Bairro Cidade Jardim, São José dos Pinhais Paraná - Brasil - CEP 83040-210 · Phone +55 41 2104-1314 · [www.sysmex.com.br](http://www.sysmex.com.br)

Program availability varies by location. Programs and specification subject to change without notice.



SAME-VISIT CBC =  
NEW POSSIBILITIES



Sysmex XW-100 Service Center  
7 am to 7 pm CT excluding holidays  
**1.800.779.7639**  
**XWServiceCenter@Sysmex.com**