

CellaVision® Proficiency Software

Digital Cell Morphology by CellaVision



How does it work?

- The software is designed to evaluate staff member's ability to perform blood and body fluid differentials
- The web-based tool can be accessed from any computer connected to the internet, using secure log-in to personalized user accounts
- The examiner can create test scenarios by uploading digital slides, either from their own CellaVision DM database or from a library of digital slides supplied by CellaVision
- The participants are automatically notified when a new test is available and instructed to perform differentials on one or more digital slides
- Upon completion, participants can browse their test result, compare their differential to that of the examiner and to those done by their co-participants
- The examiner is automatically provided with a performance report which highlights areas for future training and competency development

Why CellaVision Proficiency Software?

- Suitable for any size organization
- Can be tailored to explore identified competency gaps or meet particular training needs and goals
- Tests are quick and easy to create and manage
- Tests are quick and easy for participants to access and complete
- Results are automatically compiled
- Easy to revisit previous tests in order to keep track of participants progress

Features:

- Web-based, with instant web browser access from anywhere
- Integrated notification mechanism to participants
- Automated reporting to examiner
- Automated feedback provision to participants



Technical specification:

Recommended PC Specifications

Browser:

Latest version of Google Chrome, Microsoft Edge (Chromium based), Mozilla Firefox, Apple Safari

2 GB RAM

8 Mbps bandwidth

Screen resolution 1920x1080

Dependencies

Transfer Tool (can be used with CellaVision Remote Review Software Version 4.0 or higher)

Licences

Choose from a range of license options:

- 10, 30, 50, 70 or 100 users
- 3, 5 or 7 years

Sysmex Corporation

1-5-1-Wakinohama-Kaigandori Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0521 · www.sysmex.co.jp

Sysmex America, Inc.

577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A. · Phone +1 800 379-7639 · www.sysmex.com/us

Sysmex Canada, Inc.

5700 Explorer Drive Suite 200, Mississauga, ON L4W0C6 Canada · Phone +1 905 366-7900 · www.sysmex.ca

Sysmex Latin America and the Caribbean

Rua Joaquim Nabuco 615 - Bairro Cidade Jardim, São José dos Pinhais Paraná - Brasil - CEP 83040-210 · Phone +55 41 2104-1314 · www.sysmex.com.br