

eSupply

Inventory Management Solution





Set Your Inventory to Autopilot

Sysmex's eSupply is an innovative, automated reagent inventory management solution for smart labs.

eSupply utilizes predictive learning algorithms to proactively analyze and adjust inventory levels according to your laboratory's reagent usage patterns. It can help drive operational efficiencies across your supply chain by reducing reagent waste and limiting unnecessary shipping costs.

Smart inventory has never been so easy

eSupply automates inventory ordering processes to reduce cost, manual documentation and order frequency. With an intuitive interface and easy onboarding, eSupply is efficient, intelligent and simple to use.

SAVE

- Reduce on-hand inventory
- Consolidate and decrease shipping costs
- Lower potential costs with fewer redundant urgent or inaccurate orders

SMART

- Automatically create, notify, manage orders
- No more forgetting to place an order
- Always have the right reagent in stock

SIMPLE


- Intuitive user interface
- Easy onboarding
- Simple to manage and add users
- Accessible resources to manage process
- No installation necessary

Intelligent algorithms make it as easy as 1-2-3

1 ENTER



Technologist notices reagent inventory is running low on an analyzer



Reagent is scanned and replaced on the analyzer



eSupply automatically adds the reagent usage to the shopping list

2 REVIEW



eSupply calculates and creates your recommended order



eSupply notifies you that your order is ready to review and approve



Review, adjust if needed and approve your order

3 RECEIVE



Sysmex ships your order



Process repeats approximately every 30 days



Requirements

- Registered Customer Resource Center (CRC) user
- SNCS Edge Client: SNCS 3.5
- Google Chrome v.81 and higher
- Order reagents directly from Sysmex
- Authority to review, adjust and provide purchase order

Available for the following:

- Sysmex SP-10™, SP-50™, XE-2100™, XE-5000™, XN-10™, XN-11™, XN-20™, XN-21™, XN-350™, XN-450™, XN-550™, XS-1000i™, XT-1800i™, XT-2000i™ and XT-4000i™, RU-20™ and RU-21™, UF-5000™, UD-10™, CLINITEK Novus® Automated Urine Chemistry Analyzer

To learn more, contact your Sysmex representative

Sysmex Corporation

1-5-1-Wakinohama-Kaigandori Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0521 · www.sysmex.co.jp

Sysmex America, Inc.

577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A. · Phone +1 800 379-7639 · www.sysmex.com/us

Sysmex Canada, Inc.

5700 Explorer Drive Suite 200, Mississauga, ON L4W0C6 Canada · Phone +1 905 366-7900 · www.sysmex.ca

Sysmex Latin America and the Caribbean

Rua Joaquim Nabuco 615 - Bairro Cidade Jardim, São José dos Pinhais Paraná – Brasil – CEP 83040-210 · Phone +55 41 2104-1314 · www.sysmex.com.br