


Slide Maker by CellaVision®

HemaPrep®



High-quality smears are essential to providing accurate WBC differential results. The HemaPrep semi-automated slide maker by CellaVision simplifies, standardizes and perfects the smear making process in low volume labs. The HemaPrep is small and portable, and the simple design makes it easy to control the length and thickness of your smears.



Dimensions 90mm (W) x 155mm (H) x 135mm (D)

Order Number XU-10012

Recommended Slide Size 26 mm x 76 mm

Spreader Blade Part Number: XU-10247; 2 pcs.

Sysmex Corporation

1-5-1-Wakinohama-Kaigandori Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0521 · www.sysmex.co.jp

Sysmex America, Inc.

577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A. · Phone +1 888 879-7639 · www.sysmex.com/us

Sysmex Canada, Inc.

2680 Skymark Ave, Suite 1001, Mississauga, ON L4W 5L6, Canada · Phone +1 905 366-7900 · www.sysmex.ca

Sysmex Latin America and the Caribbean

Rua Joaquim Nabuco 615 - Bairro Cidade Jardim, São José dos Pinhais Paraná - Brasil - CEP 83040-210 · Phone +55 41 2104-1314 · www.sysmex.com.br