

A new generation of liquid QC

Chek-Stix® Urinalysis Liquid Quality Control



Quality control made effortless

Chek-Stix Liquid QC urinalysis kit is a bi-level, ready to use quality control solution for your CLINITEK Novus® Automated Urine Chemistry Analyzer. With no preparation required, laboratory professionals can spend less time on QC activities and more time running patient samples. Streamline your daily operations with Chek-Stix liquid QC.

Key Advantages



Ready to use solution



30 day open stability at room temperature



Validated for use with CLINITEK Novus® Automated Urine Chemistry Analyzer

Chek-Stix meets your lab's individual needs

Chek-Stix Liquid QC



Analyte Profile

Chemistry

pH	Urobilinogen
Protein	Specific Gravity
Glucose	Leukocyte Esterase
Ketone	
Bilirubin	
Blood (Hemoglobin)	
Nitrite	



System Compatibility

CLINITEK Novus® Automated Urine Chemistry Analyzer

Ordering Information

SMN	Each kit contains
11561695	Level 1 (Negative) Control 1 X 120 mL Level 2 (Positive) Control 1 X 120 mL

Sysmex Corporation

1-5-1-Wakino-hama-Kaigandori Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0521 · www.sysmex.co.jp

Sysmex America, Inc.

577 Aptakisic Road, Lincolnshire, IL 60069, U.S.A. · Phone +1800 379-7639 · www.sysmex.com/us

Sysmex Canada, Inc.

5700 Explorer Drive Suite 200, Mississauga, ON L4W0C6 Canada · Phone +1 905 366-7900 · www.sysmex.ca

Sysmex Latin America and the Caribbean

Rua Joaquim Nabuco 615 - Bairro Cidade Jardim, São José dos Pinhais Paraná - Brasil - CEP 83040-210 · Phone +55 41 2104-1314 · www.sysmex.com.br