



EIGHTCHECK-3WP X-TRA™

HEMATOLOGY CONTROL

XP-SERIES

Parameters	Low Abnormal			Normal			High Abnormal		
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR
	Insight Due Date #1 11/10/2019			Insight Due Date #2 12/13/2019			Insight Due Date #3 1/15/2020		
	Expiration Date: 2020-01-15			Expiration Date: 2020-01-15			Expiration Date: 2020-01-15		
	Lot Number: 92810710			Lot Number: 92810711			Lot Number: 92810712		
WBC x10 ³ /μL	3.1	0.5		6.9	0.7		17.5	1.8	
RBC x10 ⁶ /μL	2.33	0.20		4.36	0.26		5.24	0.26	
HGB g/dL	5.7	0.4		11.7	0.8		16.0	0.9	
HCT %	16.6	1.5		33.1	2.4		44.1	3.6	
MCV fL	71.2	13.0		75.9	10.5		84.2	11.3	
MCH pg	24.5	4.2		26.8	3.7		30.5	3.3	
MCHC g/dL	34.3	5.5		35.3	4.9		36.3	5.0	
PLT x10 ⁹ /μL	56	17		231	37		542	63	
Lymph (W-SCR) %	25.5	3.7		29.8	3.1		35.6	3.6	
Mixed (W-MCR) %	10.2	4.8		11.2	2.8		16.0	2.6	
Neut (W-LCR) %	64.3	5.1		59.0	3.5		48.4	3.9	
Lymph (W-SCC) x10 ³ /μL	0.8	0.2		2.1	0.5		6.2	0.8	
Mixed (W-MCC) x10 ³ /μL	0.3	0.2		0.8	0.2		2.8	0.5	
Neut (W-LCC) x10 ³ /μL	2.0	0.4		4.0	0.6		8.5	1.0	
W-SMV fL	64.9	9.6		64.3	7.0		64.5	6.4	
W-LMV fL	203.5	15.3		204.2	15.3		204.6	14.3	
RDW-SD fL	38.8	5.8		35.9	5.4		36.2	5.4	
RDW-CV %	13.4	1.7		11.3	1.6		10.1	1.7	
MPV fL	9.5	1.4		9.5	1.0		9.6	0.8	

Refer to package insert for Sysmex recommendations on setting QC targets and limits