



## XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	<b>2021-10-10</b>	Quality Control Data Due Date #1:	30-Aug-21	Quality Control Data Due Date #2:	10-Oct-21
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	12001301	12001302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 <sup>3</sup> /μL)	0.073	0.042 - 0.103	0.291	0.209 - 0.372
RBC-BF (10 <sup>6</sup> /μL)	0.024	0.020 - 0.029	0.074	0.064 - 0.085
TC-BF# (10 <sup>3</sup> /μL)	0.073	0.042 - 0.103	0.291	0.209 - 0.372
MN# (10 <sup>3</sup> /μL)	0.015	0.008 - 0.022	0.064	0.035 - 0.093
PMN# (10 <sup>3</sup> /μL)	0.057	0.026 - 0.088	0.226	0.154 - 0.299
MN%	21.2	13.5 - 28.8	22.1	13.0 - 31.2
PMN%	78.8	62.3 - 95.4	77.9	67.0 - 88.8

Please refer to the package insert for the handling of XN CHECK BF.



## XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	<b>2021-10-10</b>	Quality Control Data Due Date #1:	30-Aug-21	Quality Control Data Due Date #2:	10-Oct-21
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	12001301	12001302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 <sup>3</sup> /μL)	0.074	0.043 - 0.106	0.298	0.214 - 0.381
RBC-BF (10 <sup>6</sup> /μL)	0.024	0.020 - 0.029	0.076	0.065 - 0.087
TC-BF# (10 <sup>3</sup> /μL)	0.074	0.043 - 0.106	0.298	0.214 - 0.381
MN# (10 <sup>3</sup> /μL)	0.018	0.010 - 0.027	0.075	0.041 - 0.108
PMN# (10 <sup>3</sup> /μL)	0.056	0.026 - 0.086	0.223	0.152 - 0.295
MN%	24.7	15.8 - 33.6	25.0	14.8 - 35.3
PMN%	75.3	59.5 - 91.1	75.0	64.5 - 85.4

Please refer to the package insert for the handling of XN CHECK BF.



# XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	<b>2021-10-10</b>	Quality Control Data Due Date #1:	<b>30-Aug-21</b>	Quality Control Data Due Date #2:	<b>10-Oct-21</b>
------------------	-------------------	--------------------------------------	------------------	--------------------------------------	------------------

Lot Number:	<b>12001301</b>	<b>12001302</b>
-------------	-----------------	-----------------

Control:	<b>L1:Level 1</b>	<b>L2:Level 2</b>
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 <sup>3</sup> /μL)	0.075	0.044 - 0.107	0.301	0.217 - 0.385
RBC-BF (10 <sup>6</sup> /μL)	0.025	0.020 - 0.029	0.076	0.065 - 0.087
TC-BF# (10 <sup>3</sup> /μL)	0.075	0.044 - 0.107	0.301	0.217 - 0.385
MN# (10 <sup>3</sup> /μL)	0.018	0.010 - 0.026	0.070	0.038 - 0.101
PMN# (10 <sup>3</sup> /μL)	0.057	0.026 - 0.088	0.231	0.157 - 0.305
MN%	23.5	15.0 - 31.9	23.2	13.7 - 32.7
PMN%	76.5	60.4 - 92.6	76.8	66.1 - 87.6

Please refer to the package insert for the handling of XN CHECK BF.