



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2021-12-05	Quality Control Data Due Date #1:	25-Oct-21	Quality Control Data Due Date #2:	05-Dec-21
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	12561301	12561302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.075	0.044 - 0.107	0.290	0.209 - 0.371
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.072	0.062 - 0.083
TC-BF# (10 ³ /μL)	0.075	0.044 - 0.107	0.290	0.209 - 0.371
MN# (10 ³ /μL)	0.020	0.011 - 0.028	0.075	0.041 - 0.109
PMN# (10 ³ /μL)	0.056	0.026 - 0.086	0.215	0.146 - 0.284
MN%	26.0	16.7 - 35.4	25.8	15.2 - 36.4
PMN%	74.0	58.4 - 89.5	74.2	63.8 - 84.6

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2021-12-05	Quality Control Data Due Date #1:	25-Oct-21	Quality Control Data Due Date #2:	05-Dec-21
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	12561301	12561302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.077	0.045 - 0.110	0.295	0.212 - 0.378
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.073	0.063 - 0.084
TC-BF# (10 ³ /μL)	0.077	0.045 - 0.110	0.295	0.212 - 0.378
MN# (10 ³ /μL)	0.021	0.011 - 0.030	0.078	0.043 - 0.114
PMN# (10 ³ /μL)	0.057	0.026 - 0.087	0.217	0.147 - 0.286
MN%	26.9	17.2 - 36.5	26.6	15.7 - 37.5
PMN%	73.1	57.8 - 88.5	73.4	63.1 - 83.7

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2021-12-05	Quality Control Data Due Date #1:	25-Oct-21	Quality Control Data Due Date #2:	05-Dec-21
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	12561301	12561302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.078	0.045 - 0.110	0.297	0.214 - 0.380
RBC-BF (10 ⁶ /μL)	0.025	0.020 - 0.030	0.074	0.063 - 0.084
TC-BF# (10 ³ /μL)	0.078	0.045 - 0.110	0.297	0.214 - 0.380
MN# (10 ³ /μL)	0.020	0.011 - 0.030	0.076	0.042 - 0.110
PMN# (10 ³ /μL)	0.057	0.026 - 0.088	0.221	0.150 - 0.292
MN%	26.2	16.8 - 35.6	25.6	15.1 - 36.1
PMN%	73.8	58.3 - 89.3	74.4	64.0 - 84.8

Please refer to the package insert for the handling of XN CHECK BF.