



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2022-05-22	Quality Control Data Due Date #1:	11-Apr-22	Quality Control Data Due Date #2:	22-May-22
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	20591301	20591302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.075	0.044 - 0.107	0.294	0.211 - 0.376
RBC-BF (10 ⁶ /μL)	0.025	0.020 - 0.030	0.075	0.064 - 0.085
TC-BF# (10 ³ /μL)	0.075	0.044 - 0.107	0.294	0.211 - 0.376
MN# (10 ³ /μL)	0.020	0.011 - 0.029	0.076	0.042 - 0.110
PMN# (10 ³ /μL)	0.056	0.026 - 0.086	0.218	0.148 - 0.288
MN%	26.3	16.8 - 35.8	25.7	15.2 - 36.3
PMN%	73.7	58.2 - 89.2	74.3	63.9 - 84.7

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2022-05-22	Quality Control Data Due Date #1:	11-Apr-22	Quality Control Data Due Date #2:	22-May-22
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	20591301	20591302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.075	0.044 - 0.107	0.293	0.211 - 0.375
RBC-BF (10 ⁶ /μL)	0.026	0.021 - 0.031	0.076	0.066 - 0.087
TC-BF# (10 ³ /μL)	0.075	0.044 - 0.107	0.293	0.211 - 0.375
MN# (10 ³ /μL)	0.020	0.011 - 0.030	0.080	0.044 - 0.116
PMN# (10 ³ /μL)	0.055	0.025 - 0.085	0.213	0.145 - 0.281
MN%	27.1	17.4 - 36.9	27.4	16.2 - 38.6
PMN%	72.9	57.6 - 88.2	72.6	62.4 - 82.8

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2022-05-22	Quality Control Data Due Date #1:	11-Apr-22	Quality Control Data Due Date #2:	22-May-22
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	20591301	20591302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.076	0.044 - 0.108	0.299	0.215 - 0.382
RBC-BF (10 ⁶ /μL)	0.026	0.021 - 0.030	0.077	0.066 - 0.088
TC-BF# (10 ³ /μL)	0.076	0.044 - 0.108	0.299	0.215 - 0.382
MN# (10 ³ /μL)	0.020	0.011 - 0.029	0.079	0.043 - 0.114
PMN# (10 ³ /μL)	0.056	0.026 - 0.086	0.220	0.150 - 0.290
MN%	26.6	17.0 - 36.1	26.4	15.6 - 37.2
PMN%	73.4	58.0 - 88.8	73.6	63.3 - 83.9

Please refer to the package insert for the handling of XN CHECK BF.