



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2022-07-17	Quality Control Data Due Date #1:	06-Jun-22	Quality Control Data Due Date #2:	17-Jul-22
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	21151301	21151302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.076	0.044 - 0.108	0.324	0.233 - 0.415
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.072	0.062 - 0.083
TC-BF# (10 ³ /μL)	0.076	0.044 - 0.108	0.324	0.233 - 0.415
MN# (10 ³ /μL)	0.020	0.011 - 0.029	0.084	0.046 - 0.122
PMN# (10 ³ /μL)	0.056	0.026 - 0.086	0.240	0.163 - 0.316
MN%	26.4	16.9 - 35.8	26.0	15.3 - 36.6
PMN%	73.6	58.2 - 89.1	74.0	63.6 - 84.4

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2022-07-17	Quality Control Data Due Date #1:	06-Jun-22	Quality Control Data Due Date #2:	17-Jul-22
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	21151301	21151302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.076	0.044 - 0.108	0.325	0.234 - 0.416
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.073	0.063 - 0.083
TC-BF# (10 ³ /μL)	0.076	0.044 - 0.108	0.325	0.234 - 0.416
MN# (10 ³ /μL)	0.021	0.011 - 0.030	0.088	0.049 - 0.128
PMN# (10 ³ /μL)	0.055	0.026 - 0.085	0.237	0.161 - 0.313
MN%	27.1	17.3 - 36.8	27.1	16.0 - 38.3
PMN%	72.9	57.6 - 88.2	72.9	62.7 - 83.1

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2022-07-17	Quality Control Data Due Date #1:	06-Jun-22	Quality Control Data Due Date #2:	17-Jul-22
------------------	-------------------	--------------------------------------	------------------	--------------------------------------	------------------

Lot Number:	21151301	21151302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.078	0.045 - 0.110	0.333	0.240 - 0.426
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.074	0.064 - 0.084
TC-BF# (10 ³ /μL)	0.078	0.045 - 0.110	0.333	0.240 - 0.426
MN# (10 ³ /μL)	0.020	0.011 - 0.029	0.085	0.047 - 0.124
PMN# (10 ³ /μL)	0.058	0.026 - 0.089	0.248	0.168 - 0.327
MN%	25.8	16.5 - 35.1	25.6	15.1 - 36.2
PMN%	74.2	58.6 - 89.8	74.4	63.9 - 84.8

Please refer to the package insert for the handling of XN CHECK BF.