



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2020-07-19	Quality Control Data Due Date #1:	08-Jun-20	Quality Control Data Due Date #2:	19-Jul-20
------------------	-------------------	--------------------------------------	------------------	--------------------------------------	------------------

Lot Number:	01181301	01181302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.075	0.044 - 0.107	0.312	0.225 - 0.400
RBC-BF (10 ⁶ /μL)	0.023	0.019 - 0.028	0.073	0.063 - 0.083
TC-BF# (10 ³ /μL)	0.075	0.044 - 0.107	0.312	0.225 - 0.400
MN# (10 ³ /μL)	0.018	0.010 - 0.025	0.072	0.040 - 0.105
PMN# (10 ³ /μL)	0.058	0.027 - 0.089	0.240	0.163 - 0.317
MN%	23.3	14.9 - 31.7	23.2	13.7 - 32.7
PMN%	76.7	60.6 - 92.8	76.8	66.1 - 87.6

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2020-07-19	Quality Control Data Due Date #1:	08-Jun-20	Quality Control Data Due Date #2:	19-Jul-20
------------------	-------------------	--------------------------------------	------------------	--------------------------------------	------------------

Lot Number:	01181301	01181302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.077	0.045 - 0.110	0.317	0.228 - 0.406
RBC-BF (10 ⁶ /μL)	0.024	0.019 - 0.029	0.074	0.064 - 0.084
TC-BF# (10 ³ /μL)	0.077	0.045 - 0.110	0.317	0.228 - 0.406
MN# (10 ³ /μL)	0.021	0.011 - 0.030	0.087	0.048 - 0.126
PMN# (10 ³ /μL)	0.056	0.026 - 0.087	0.230	0.156 - 0.303
MN%	26.9	17.2 - 36.6	27.5	16.2 - 38.8
PMN%	73.1	57.8 - 88.5	72.5	62.4 - 82.7

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2020-07-19	Quality Control Data Due Date #1:	08-Jun-20	Quality Control Data Due Date #2:	19-Jul-20
------------------	-------------------	--------------------------------------	------------------	--------------------------------------	------------------

Lot Number:	01181301	01181302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.079	0.046 - 0.112	0.326	0.234 - 0.417
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.074	0.063 - 0.084
TC-BF# (10 ³ /μL)	0.079	0.046 - 0.112	0.326	0.234 - 0.417
MN# (10 ³ /μL)	0.019	0.011 - 0.028	0.080	0.044 - 0.115
PMN# (10 ³ /μL)	0.059	0.027 - 0.091	0.246	0.167 - 0.325
MN%	24.7	15.8 - 33.6	24.4	14.4 - 34.4
PMN%	75.3	59.5 - 91.1	75.6	65.0 - 86.2

Please refer to the package insert for the handling of XN CHECK BF.