



XN CHECK BF™

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2019-10-13	Quality Control Data Due Date #1:	2-Sep-19	Quality Control Data Due Date #2:	13-Oct-19
------------------	-------------------	--------------------------------------	----------	--------------------------------------	-----------

Lot Number:	92031301	92031302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF($10^3/\mu\text{L}$)	0.078	0.045 - 0.111	0.310	0.223 - 0.397
RBC-BF($10^6/\mu\text{L}$)	0.026	0.021 - 0.031	0.076	0.065 - 0.087
TC-BF#($10^3/\mu\text{L}$)	0.078	0.045 - 0.111	0.310	0.223 - 0.397
MN#($10^3/\mu\text{L}$)	0.020	0.011 - 0.029	0.078	0.043 - 0.113
PMN#($10^3/\mu\text{L}$)	0.059	0.027 - 0.091	0.233	0.158 - 0.308
MN%	25.0	16.0 - 34.0	25.0	14.8 - 35.3
PMN%	75.0	59.3 - 90.8	75.0	64.5 - 85.5

Please refer to the package insert for the handling of XN CHECK BF.