



XN CHECK BF™

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2019-10-13	Quality Control Data Due Date #1:	2-Sep-19	Quality Control Data Due Date #2:	13-Oct-19
------------------	-------------------	--------------------------------------	----------	--------------------------------------	-----------

Lot Number:	92031301	92031302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF($10^3/\mu\text{L}$)	0.078	0.045 - 0.111	0.313	0.225 - 0.401
RBC-BF($10^6/\mu\text{L}$)	0.026	0.021 - 0.031	0.076	0.065 - 0.087
TC-BF#($10^3/\mu\text{L}$)	0.078	0.045 - 0.111	0.313	0.225 - 0.401
MN#($10^3/\mu\text{L}$)	0.020	0.011 - 0.029	0.082	0.045 - 0.119
PMN#($10^3/\mu\text{L}$)	0.058	0.027 - 0.089	0.231	0.157 - 0.305
MN%	26.0	16.6 - 35.4	26.3	15.5 - 37.1
PMN%	74.0	58.5 - 89.5	73.7	63.4 - 84.0

Please refer to the package insert for the handling of XN CHECK BF.