



EIGHTCHECK-3WP X-TRA™

HEMATOLOGY CONTROL

XP-SERIES

Low Abnormal	Normal	High Abnormal
<i>Insight</i> Due Date #1: 11/10/2019	<i>Insight</i> Due Date #2: 12/13/2019	<i>Insight</i> Due Date #3: 1/15/2020
Expiration Date: 2020-01-15	Expiration Date: 2020-01-15	Expiration Date: 2020-01-15
Lot Number: 92820710	Lot Number: 92820711	Lot Number: 92820712

Parameters	Low Abnormal			Normal			High Abnormal		
	Mean	Limit	BCR	Mean	Limit	BCR	Mean	Limit	BCR
WBC x10 ³ /μL	3.2	0.5		6.9	0.7		17.5	1.8	
RBC x10 ⁶ /μL	2.34	0.20		4.35	0.26		5.25	0.26	
HGB g/dL	5.7	0.4		12.1	0.8		15.7	0.9	
HCT %	16.6	1.5		34.0	2.4		43.4	3.6	
MCV fL	70.9	13.0		78.2	10.5		82.7	11.3	
MCH pg	24.4	4.2		27.8	3.7		29.9	3.3	
MCHC g/dL	34.3	5.5		35.6	4.9		36.2	5.0	
PLT x10 ³ /μL	58	17		234	37		533	63	
Lymph (W-SCR) %	25.9	3.7		29.0	3.1		34.9	3.6	
Mixed (W-MCR) %	10.9	4.8		11.5	2.8		16.2	2.6	
Neut (W-LCR) %	63.2	5.1		59.5	3.5		48.9	3.9	
Lymph (W-SCC) x10 ³ /μL	0.8	0.2		2.0	0.5		6.1	0.8	
Mixed (W-MCC) x10 ³ /μL	0.3	0.2		0.8	0.2		2.8	0.5	
Neut (W-LCC) x10 ³ /μL	2.1	0.4		4.1	0.6		8.6	1.0	
W-SMV fL	64.5	9.6		64.4	7.0		64.7	6.4	
W-LMV fL	204.3	15.3		204.9	15.3		204.9	14.3	
RDW-SD fL	36.5	5.8		33.7	5.4		35.0	5.4	
RDW-CV %	12.4	1.7		9.9	1.6		10.0	1.7	
MPV fL	9.9	1.4		9.8	1.0		9.9	0.8	

Refer to package insert for Sysmex recommendations on setting QC targets and limits