



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2020-05-24	Quality Control Data Due Date #1:	13-Apr-20	Quality Control Data Due Date #2:	24-May-20
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	00621301	00621302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.077	0.045 - 0.109	0.308	0.222 - 0.394
RBC-BF (10 ⁶ /μL)	0.026	0.021 - 0.031	0.076	0.065 - 0.087
TC-BF# (10 ³ /μL)	0.077	0.045 - 0.109	0.308	0.222 - 0.394
MN# (10 ³ /μL)	0.018	0.010 - 0.026	0.072	0.040 - 0.104
PMN# (10 ³ /μL)	0.059	0.027 - 0.091	0.236	0.160 - 0.312
MN%	22.9	14.7 - 31.1	23.4	13.8 - 33.0
PMN%	77.1	60.9 - 93.3	76.6	65.9 - 87.3

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2020-05-24	Quality Control Data Due Date #1:	13-Apr-20	Quality Control Data Due Date #2:	24-May-20
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	00621301	00621302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.075	0.044 - 0.107	0.315	0.227 - 0.403
RBC-BF (10 ⁶ /μL)	0.026	0.021 - 0.031	0.077	0.066 - 0.088
TC-BF# (10 ³ /μL)	0.075	0.044 - 0.107	0.315	0.227 - 0.403
MN# (10 ³ /μL)	0.020	0.011 - 0.029	0.086	0.047 - 0.125
PMN# (10 ³ /μL)	0.055	0.025 - 0.085	0.229	0.156 - 0.302
MN%	26.3	16.8 - 35.8	27.4	16.2 - 38.6
PMN%	73.7	58.2 - 89.2	72.6	62.4 - 82.8

Please refer to the package insert for the handling of XN CHECK BF.



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2020-05-24	Quality Control Data Due Date #1:	13-Apr-20	Quality Control Data Due Date #2:	24-May-20
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	00621301	00621302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.079	0.046 - 0.112	0.317	0.228 - 0.406
RBC-BF (10 ⁶ /μL)	0.026	0.021 - 0.031	0.076	0.065 - 0.087
TC-BF# (10 ³ /μL)	0.079	0.046 - 0.112	0.317	0.228 - 0.406
MN# (10 ³ /μL)	0.020	0.011 - 0.029	0.080	0.044 - 0.116
PMN# (10 ³ /μL)	0.059	0.027 - 0.091	0.237	0.161 - 0.313
MN%	25.7	16.4 - 35.0	25.1	14.8 - 35.4
PMN%	74.3	58.7 - 89.9	74.9	64.4 - 85.4

Please refer to the package insert for the handling of XN CHECK BF.