



XN CHECK BF™

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

| | | | | | |
|------------------|-------------------|--------------------------------------|----------|--------------------------------------|-----------|
| Expiration Date: | 2019-08-18 | Quality Control Data Due Date #1: | 8-Jul-19 | Quality Control Data Due Date #2: | 18-Aug-19 |
|------------------|-------------------|--------------------------------------|----------|--------------------------------------|-----------|

| | | |
|-------------|-----------------|-----------------|
| Lot Number: | 91471301 | 91471302 |
|-------------|-----------------|-----------------|

| | | |
|----------|-------------------|-------------------|
| Control: | L1:Level 1 | L2:Level 2 |
|----------|-------------------|-------------------|

| PARAMETERS | MEAN | EXPECTED RANGE | MEAN | EXPECTED RANGE |
|------------------------------|-------|----------------|-------|----------------|
| WBC-BF($10^3/\mu\text{L}$) | 0.082 | 0.048 - 0.116 | 0.333 | 0.240 - 0.426 |
| RBC-BF($10^6/\mu\text{L}$) | 0.026 | 0.021 - 0.031 | 0.076 | 0.065 - 0.087 |
| TC-BF#($10^3/\mu\text{L}$) | 0.082 | 0.048 - 0.116 | 0.333 | 0.240 - 0.426 |
| MN#($10^3/\mu\text{L}$) | 0.020 | 0.011 - 0.029 | 0.079 | 0.052 - 0.106 |
| PMN#($10^3/\mu\text{L}$) | 0.062 | 0.029 - 0.095 | 0.254 | 0.173 - 0.335 |
| MN% | 23.9 | 15.3 - 32.5 | 23.6 | 18.6 - 28.6 |
| PMN% | 76.1 | 60.1 - 92.1 | 76.4 | 65.7 - 87.1 |

Please refer to the package insert for the handling of XN CHECK BF.