



# XN CHECK BF™

## BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

<b>Expiration Date:</b>	<b>2019-08-18</b>	<b>Quality Control Data Due Date #1:</b>	<b>8-Jul-19</b>	<b>Quality Control Data Due Date #2:</b>	<b>18-Aug-19</b>
-------------------------	-------------------	--	-----------------	--	------------------

<b>Lot Number:</b>	<b>91471301</b>	<b>91471302</b>
--------------------	-----------------	-----------------

<b>Control:</b>	<b>L1:Level 1</b>	<b>L2:Level 2</b>
-----------------	-------------------	-------------------

<b>PARAMETERS</b>	<b>MEAN</b>	<b>EXPECTED RANGE</b>	<b>MEAN</b>	<b>EXPECTED RANGE</b>
WBC-BF( $10^3/\mu\text{L}$ )	0.079	0.046 - 0.112	0.326	0.235 - 0.417
RBC-BF( $10^6/\mu\text{L}$ )	0.026	0.021 - 0.031	0.076	0.065 - 0.087
TC-BF#( $10^3/\mu\text{L}$ )	0.079	0.046 - 0.112	0.326	0.235 - 0.417
MN#( $10^3/\mu\text{L}$ )	0.020	0.011 - 0.029	0.082	0.054 - 0.110
PMN#( $10^3/\mu\text{L}$ )	0.059	0.027 - 0.091	0.244	0.166 - 0.322
MN%	25.6	16.4 - 34.8	25.2	19.9 - 30.5
PMN%	74.4	58.8 - 90.0	74.8	64.3 - 85.3

Please refer to the package insert for the handling of XN CHECK BF.