



XN CHECK BF™

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2019-10-13	Quality Control Data Due Date #1:	2-Sep-19	Quality Control Data Due Date #2:	13-Oct-19
------------------	-------------------	--------------------------------------	----------	--------------------------------------	-----------

Lot Number:	92031301	92031302
-------------	-----------------	-----------------

Control:	L1:Level 1	L2:Level 2
----------	-------------------	-------------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF($10^3/\mu\text{L}$)	0.081	0.047 - 0.115	0.321	0.231 - 0.411
RBC-BF($10^6/\mu\text{L}$)	0.026	0.021 - 0.031	0.076	0.065 - 0.087
TC-BF#($10^3/\mu\text{L}$)	0.081	0.047 - 0.115	0.321	0.231 - 0.411
MN#($10^3/\mu\text{L}$)	0.021	0.012 - 0.030	0.083	0.046 - 0.120
PMN#($10^3/\mu\text{L}$)	0.060	0.028 - 0.092	0.238	0.162 - 0.314
MN%	25.7	16.4 - 35.0	25.9	15.3 - 36.5
PMN%	74.3	58.7 - 89.9	74.1	63.7 - 84.5

Please refer to the package insert for the handling of XN CHECK BF.