




## BEYOND A BETTER BOX™

### DI-60 REMOTE REVIEW SOFTWARE

*Review Slide Files and Cell Images Anytime, Anywhere*



# DI-60 REMOTE REVIEW SOFTWARE

By using the DI-60 Remote Review Software, clinicians obtain real-time access to analyzed slides from their office, clinic or home. In addition to today's cell images, access is available to a patient's historical images for comparison.

## LABORATORY



When used on the laboratory's network, the software can significantly enhance workflow and labor resources. Multiple sites can now centralize morphology expertise. Additionally, the organization can achieve staffing flexibility since it is now possible for any lab on a network to review files and images.



## ALTERNATE SITE




## OFFICE

# DI-60 REMOTE REVIEW SOFTWARE

## Product Information

DI-60 Remote Review Software provides remote users access to analyzed slides and the ability to reclassify cells and sign off slides from another location. This allows for multiple work-stations within the lab, as well as real-time consultation by a supervisor or clinician from within or outside the hospital.



## KEY FEATURES

- Verify, correct or collaborate on analyzed blood, body fluids and digital slides remotely.
- Access patient's cell images, both current and historical, remotely.
- Identify and comment on any number of tagged Regions of Interest when reviewing body fluids or other specimens.
- Collaborate within/between hospitals in real time.
- View cell counter results and flags.

## Technical Specifications

### Recommended PC Specifications

- Operating systems: Win Vista/Win 7
- 256 MB graphics RAM
- Ethernet adapter 10/100 Mbps
- 512 MB RAM
- 1 GB free disk space
- CPU Pentium IV or higher

### Supported Instrument

- Sysmex DI-60

### Supported Applications

- DI-60 Peripheral Blood Application
- DI-60 Body Fluid Application

### DI-60 Remote Review Software

- Single Workstation on LAN
- Part number: BC188562



---

## BEYOND A BETTER BOX™

Sysmex delivers a comprehensive and intuitive ecosystem that improves performance with next generation diagnostics, advanced tools and technologies, process optimization and BeyondCare<sup>SM</sup> harmonized support.

---

**Sysmex Corporation**  
1-5-1-Wakinohama-Kaigandori  
Chuo-ku, Kobe 651-0073, Japan  
Tel. +81 (78) 265-0521  
[www.sysmex.co.jp](http://www.sysmex.co.jp)

**Sysmex America, Inc.**  
577 Aptakistic Road  
Lincolnshire, IL 60069, U.S.A.  
Tel. +1 (800) 379-7639  
[www.sysmex.com/us](http://www.sysmex.com/us)

**Sysmex Canada, Inc.**  
5700 Explorer Drive  
Suite 200,  
Mississauga, ON L4W0C6  
Canada  
Tel. +1 (905) 366-7900  
[www.sysmex.ca](http://www.sysmex.ca)

**Sysmex Latin America and  
the Caribbean**  
Rua Joaquim Nabuco 615 -  
Bairro Cidade Jardim  
São José dos Pinhais  
Paraná – Brasil – CEP 83040-210  
Tel. +55 (41) 2104-1314  
[www.sysmex.com/la](http://www.sysmex.com/la)

